

## Lucent AD (p-Tau 181), K2 EDTA Plasma

Test ID: 104142 p-Tau 181

Notification Date: 06Jul2023

Effective Date: 06Jul2023

Explanation: On the effective date, the p-Tau 181 test, performed by the Quanterix Corporation Lab, will be updated with the following changes.

Current Test Name
p-Tau 181

New Test Name
Lucent AD p-Tau 181

Current Cut-Off
$\leq 16.5$ pg/mL

New Cut-Off
$\leq 14.2$ pg/mL

Current Specimen Stability
4 days at 2-8°C 4 days at -20°C

New Specimen Stability
7 days at 2-8°C 7 days at -20°C

Questions Contact Clinical Services email: [support@lucentdiagnostics.com](mailto:support@lucentdiagnostics.com)